

CURIOX'S LAMINAR WASH TECHNOLOGY




Fig 1 Curiox's centrifuge-less cell preparation platform is enabled by our **wall-less plate** and **laminar flow washer**.

(a) MINI 1000 is a personal bench-top instrument that processes 8-wells simultaneously, on a strip plate consisting 16-wells.

(b) HT1000 washes 96 wells in a fully automated process.

(c) The Laminar Wash™ (LW) plate (96 well) or strip (16 well) consists of an array of hydrophilic spots surrounded by hydrophobic surface, which functions as a virtual wall that separates each spot. Each spot can process from a single cell to as many as 10 million cells without the mechanical stress and cell losses associated with centrifugation.

FULL AUTOMATION for CELLULAR STAINING and WHOLE BLOOD PROCESSING for FLOW CYTOMETRY ASSAYS



The **Laminar Wash™ AUTO 1000** software interface allows modifications to sample volume, buffer volume, antibody volume, number of washes, incubation time and temperature. In addition, it comes with pre-programmed protocols to ensure user consistency.

User-friendly Interface: The software prompts the user to enter the parameters of existing SOPs or protocols with no coding or programming of the liquid handling system.

- Automated** – reduces manual pipetting errors and errors associated with multiple personnel changes
- Turnkey surface and intracellular staining** – automation improves lab efficiency and reduces waste by reducing repeat experiments
- Enables compliance** – every step of the protocol is recorded in the software
- Rapid time to results** – AUTO 1000 processes 2 plates in under 2 hours using pre-programmed intracellular staining protocols
- Consistent results** – across users and locations
- Cleaner data** – Laminar Wash™ reduces the amount of background and debris, especially with tissue samples

AUTOMATED LAMINAR WASH™ PRODUCED CONSISTENT DATA with LOW COEFFICIENT of VARIANCE




Fig 2 By using Laminar Wash™ in a single wash step during immunophenotyping of human PBMCs, the user achieved **(a)** consistency and **(b)** lower coefficient of variance even between different operators. *Webinar available from www.curiox.com*

LAMINAR WASH™ ENABLED ACCURATE IDENTIFICATION of TUMOR INFILTRATING LYMPHOCYTE (TIL) POPULATION




Fig 3(a) Centrifuge wash method – cell loss and mechanical stress through sequential pelleting and resuspension **Fig 3(b)** Laminar Wash™ method – less tumor debris, but higher retention of TILs and better resolution of populations

LAMINAR WASH™ IMPROVED RETENTION and RETAINED ENDOGENOUS PROTEIN EXPRESSION




Fig 4(a) In processing of low cell numbers (50-100), LW retained higher cell count with better cellular morphology and brighter endogenous fluorescent expression. The dim protein expression is evidence of mechanic stresses caused by centrifugation. *Data from a biotech company in San Diego*

Fig 4(b) In subcutaneous murine melanoma tumor samples, LW retained 40% more live cells than centrifuge, with less tumor debris and better antibody staining.

LAMINAR WASH™ DECREASED BACKGROUND in B CELL CROSSMATCH without PRONASE TREATMENT




Fig 5(a) Decreased background in pronase-independent B cell crossmatching of PBL fresh cells with Laminar Wash™ **5(b)** ΔMCF was eliminated with Laminar Wash™ in pronase-untreated B cells bound to the wells of negative control serum. *Webinar available from www.curiox.com*

LAMINAR WASH™ LED to IMPROVED CELL VIABILITY of TILs and LESS CYTOMETER CLOGGING




Fig 6 FACS plots by **(a)** centrifuge vs. **(b)** Laminar Wash™ methods. The dissociated tumor samples from Raji-inoculated Balb/c mice showed higher viability and increased debris removal with Laminar Wash™, leading to less clogging and more consistent reading.

ACKNOWLEDGEMENTS

We would like to give special thanks to the following scientists: Dr. Arnaud Colantonio at Adicet Bio (Fig 2), Dr. Christoph Eberle at Charles River Laboratories in Worcester, MA (Fig 3), and Dr. Thomas Tan, Jonathan Hsu, and Benjamin Keller at Elstar Therapeutics, Cambridge, MA (Fig 6) for sharing their valuable insights and data.

Come see **Laminar Wash™ AUTO 1000** and your personal-sized **MINI1000** in action at Booth #207!

Learn more about how **Laminar Wash™** can benefit your **whole blood processing** assays at Poster #113.

